

Installing the ActiveX MeadCo ScriptX Print Control in Internet Explorer



The ActiveX MeadCo ScriptX Print Control is required to properly view and print lab reports.


Requirements:

Windows based PC or laptop (Windows 7, Windows Vista, Windows XP).

Internet Explorer version 6 or higher.

*Screenshots may vary depending on version of Windows and/or Internet Explorer.

If the MeadCo ScriptX Print Control is not installed the following message will appear upon login.



Page 1 of 3




Installing the ActiveX MeadCo ScriptX Print Control in Internet Explorer



This window will pop up in the background and request to install the ActiveX control.

- 1. Click on the bar near the top of the web page,
- 2. Click "Install this add on for all users of this computer."




Page 2 of 3 Rev2




Installing the ActiveX MeadCo ScriptX Print Control in Internet Explorer



You will see the circle on the Internet Explorer tab "thinking".



After ActiveX is done installing the thinking circle will turn to the Internet Explorer E. This pop-up window can now be closed.



Page 3 of 3 Rev2