







HOW TO USE BD™ Universal Viral Transport



COLLECTION AND TRANSPORT GUIDE



220527

BD™ Universal Viral Transport

3mL vial with 1 sterile nylon flocked regular and 1 sterile nylon flocked flexible minitip swabs

For the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasmas or ureaplasmas.



220531

BD™ Universal Viral Transport

3mL vial and 1 sterile nylon flocked flexible mini-tip flocked swab

For the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasmas or ureaplasmas.

See product package insert for instructions for use and manufacturer's recommendations.

Eswab is a trademark of Copan Italia S.P.A.
©2016 BD. BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.